# SUSPENSION BOOT





#### Description

Help support to relieve pressure in the heel area. Especially indicated for bedridden patients.

### Effects abd advantages

Possibility of walking short distances. The material from which it is made promotes air circulation. Support attachment via additional hook-and-loop closure at forefoot. Control of foot position and movement through the semi-rigid splint. Suitable for left and right foot.

### Composition

Blue nylon and black velor fabric. Filled with foam. Polypropylene tapes. Velcro.

### Dimensions

Universal size.

#### Maintenance

Remove the plastic ferrule. Hand or machine wash with a delicate program at a maximum temperature of 40°C. Do not use abrasive products, whiteners or bleach. Do not iron.

#### Warranty


The warranty period is 2 years. This begins



with the date of purchase from the dealer. The guarantee does not cover normal wear and tearof the product and does not replace legal guarantees. Please contact the company's quality department in case of an incident to resolve it. Product lifespan 1 year.

# Anti-decubitus





## Colors





MODEL	SIZE	REF.
Suspension boot	Universal	501310





CE

catálogo ubiotex www.ubiotex.com