

SUAPEL ELBOW PAD



Description

Sanitized suapel elbow pad.

Effects and advantages

Prevents redness and bedsores by avoiding direct contact with other surfaces.

Composition

Suapel: 60% polyester and 40% acrylic.
Buckle fastening.

Dimensions

22x34x2,5cm.

Maintenance

Hand or machine washable with neutral detergent at a maximum temperature of 40°C.
Tumble dry at reduced temperature.
Do not iron.
Do not bleach.
Dry horizontally.


Warranty

The warranty period is 2 years. It starts from the date of purchase from the dealer. The warranty does not cover normal wear and tear of the product and does not replace statutory warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Anti-decubitus

Suapel elbow pad

Models



Colours



Navy
blue




Grey



White

MODEL	DIMENSIONS (CM)	COLOUR	REF.
Suapel elbow pad	22x34x2,5	White	112100
		Navy blue	112102
		Grey	112104



antiol or



anti bacterias



anti hon gos



anti ácar os



anti CHINCHES



SUAVE



TRANSPIRABLE



LAVABLE

CE