

PADDED SUAPEL ELBOW PAD



Description

Padded suapel elbow pad.

Effects and advantages

Prevents the generation of redness and bedsores by avoiding direct contact with other surfaces. Recommended in cases where there is a risk of bedsores as, in addition to being made of suapel fabric, it has an inner filling of silicone fibre that considerably relieves pressure. In addition to preventing bedsores, it favours their treatment by eliminating pressure in the area.

Composition

- Suapel: 60% polyester and 40% acrylic.
 - Terlenka: 65% polyester and 35% cotton.
- Filled with fibre and with velcro fastening.

Dimensions

22x34x2,5cm.

Maintenance

Hand or machine washable with neutral detergent at a maximum temperature of 40°C. Tumble dry at reduced temperature. Do not iron. Do not bleach. Dry horizontally.

Warranty


The warranty period is 2 years. It starts from




the date of purchase from the dealer. The warranty does not cover normal wear and tear of the product and does not replace statutory warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Anti-decubitus

Padded suapel elbow pad Models



Colours



MODEL	DIMENSIONS (CM)	COLOUR	REF.
Padded suapel elbow pad	22x34x2,5	White	112120
		Navy blue	112201



antiol or



anti bacterias



anti hon gos



anti ácar os



anti CHINCHES



SUAVE



TRANSPIRABLE



LAVABLE

