# FIBRE CORE BACKREST AND SEAT





#### Description

Backrest with seat made of independent antidecubitus fibre cores.

# Effects and advantages

Prevents the development of bedsores in people who have to spend long periods of time in a seated position.

#### Composition

Fibrefill. Saniluxe polyurethane: 60% polyester, 40% polyurethane. Terlenka.

Elastic strap with buckle to adjust to the backrest. Straps to adjust to the seat.

# Dimensions

Seat: 44x44x16cm. Backrest: 44x54x16cm.

# Maintenance


Hand wash or machine wash with neutral detergent at 95°C. Do not spin-dry. Avoid abrasive products. Do not use bleach. Tumble dry at low temperature. Do not iron. Dry flat.

### Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

# Anti-decubitus

Fibre core backrest and seat. Models



Grosor: 16cm.

# Colours



MODEL	DIMENSIONS (CM)	REF.
Fibre core backrest and seat.	44x44x16 44x56x16	114152



