

LEG TRANSFER BELT



Description

Leg transfer belt.

Effects and benefits

Facilitates the mobilisation and transfer of the patient's lower limbs. It offers easy handling by caregivers.

Composition

Quilted fabric (20% suede, 20% velour and 60% foam).

Polypropylene webbing.

Velcro safety fastening.

Dimensions

To be consulted.

Maintenance

Hand or machine wash with a programme for delicate garments at a maximum temperature of 40°C.

Do not use abrasive or bleaching products.

Do not iron.

Do not use bleach.


Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the

company's quality department in the event of an incident in order to resolve it.

Restraints

Leg transfer belt Models



Colours



Navy
blue

MODEL	COLOUR	DIMENSIONS (CM)	REF.
Leg transfer belt	Navy blue	To be consulted	306604

