

## TRANSFER BLANKET



### Description

Mobilising and transfer blanket.

### Effects and benefits

Specially designed for patients with psychomotor difficulties. Helps to incorporate and/or transfer the patient.

### Composition

Polyurethane, polypropylene straps and neoprene handles.

### Dimensions

To be consulted

### Maintenance

Regular hand or machine washing at a maximum temperature of 40°C is sufficient.

Dry in an upright position without direct heat source.

Do not bleach.


### Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.



# Restraints

## Transfer blanket Models



MODEL	COLOUR	DIMENSIONS (CM)	REF.
Transfer blanket	Grafite	To be consulted	306600

