Restraints

STANDARD CRANE STRAP





Description Foot strap for standing lift.

Effects and benefits

Used for added safety on the footplate of the sit-to-stand lift.

Composition

Quilted fabric (20% suede, 20% velour and 60% foam). Polypropylene webbing. Velcro safety fastening.

Dimensions

Universal.

Maintenance

Hand or machine wash with a programme for delicate garments at a maximum temperature of 40°C.

Do not use abrasive or bleaching products. Do not iron. Do not use bleach.


Warranty

The warranty period is two years from the date of purchase. This warranty does not cover

normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Restraints

Standard crane strap Models



Colours



MODEL	COLOUR	DIMENSIONS (CM)	REF.
Standard crane strap	Navy blue	Universal	306710

