

LEG REST BAND



Description

Legrest band for chair or armchair.

Effects and benefits

The calf rest provides firm and comfortable support in the calf area and has a secure fastening at the back.

Indications for use

To position the article correctly, the patient should be placed comfortably in the chair or armchair and the leg support band should be passed behind the calves, adjusting it to the most comfortable height for the patient.

Composition

- Padding: 20% anteline (100% polyester), 20% velour (100% polyamide) and 60% foam.
Straps: 100% polypropylene.

Dimensiones

Universal.

Maintenance


Regular hand or machine washing at a maximum temperature of 40°C is sufficient.
Dry in a horizontal position without exposing to direct heat sources.

Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Restraints

Leg rest band Models



Colours



Navy blue

MODEL	FABRIC	CLOSURE	REF.
Leg rest band	Padded	Buckles	307201

