

## ABDOMINAL RESTRAINTS BELTS BED 105



### Description

Abdominal belt for 105 bed with buckle, magnet or ironclip fasteners.

### Effects and benefits

Total immobilisation of the waist. Prevents falls with partial immobilisation of the body focused on the abdomen area. Lateralisation of the user is permitted and, in the event that the user wishes to be immobilised face up, two polypropylene straps are available for fastening.

### Composition

- Quilted fabric: 20% suede (100% polyester), 20% velour (100% polyamide) and 60% foam.
- Canvas: 75% cotton, 25% polyester.
- Polypropylene: 100% polypropylene.

### Dimensions

S, M and L.

### Maintenance

They do not require special care. It is enough to wash them periodically by hand or machine, at a maximum temperature of 40°C.

Dry in a horizontal position without exposing to direct heat sources such as cookers or heaters.

Do not tumble dry.

Do not iron.




### Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

# Restraints

## Abdominal restraints belts bed 105

### Models



MODEL	FABRIC	CLOSURE	SIZE	REF.
Abdominal restraints belts bed 105	Polypropylene	Buckles	M	401111
		Ironclip	M	401130
			L	401132
		Magnet	M	401105
	L		401107	
	Padded	Buckles	M	401205
			L	401207
		Ironclip	M	401270
			L	401272
		Magnet	M	401245
			L	401247
	Canvas	Ironclip	M	401455
			L	401457
		Magnet	M	401410
L			401412	

### Colours



Black



Canvas

