

ABDOMINAL RESTRAINTS BELTS BED 90



Description

Abdominal belt for 90° bed with buckle, magnet or ironclip fasteners.

Effects and benefits

Total immobilisation of the waist. Prevents falls with partial immobilisation of the body focused on the abdomen area. Lateralisation of the user is permitted and, in the event that the user wishes to be immobilised face up, two polypropylene straps are available for fastening.

Composition

- Quilted fabric: 20% suede (100% polyester), 20% velour (100% polyamide) and 60% foam.
- Canvas: 75% cotton, 25% polyester.
- Polypropylene: 100% polypropylene.

Dimensions

S, M and L.

Maintenance

They do not require special care. It is enough to wash them periodically by hand or machine, at a maximum temperature of 40°C.

Dry in a horizontal position without exposing to direct heat sources such as cookers or heaters.

Do not tumble dry.

Do not iron.




Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Restraints

Abdominal restraints belts bed 90

Models



| MODEL | FABRIC | CLOSURE | SIZE | REF. |
|-----------------------------------|---------------|--------------|---------|--------|
| Abdominal restraints belts bed 90 | Polypropylene | Buckles | M | 401100 |
| | | | L | 401102 |
| | | Bucklesbasic | M | 401115 |
| | | | Magnet | M |
| | | L | | 401122 |
| | | Ironclip | M | 401150 |
| | | | L | 401152 |
| | | Padded | Buckles | S |
| | M | | | 401200 |
| | L | | | 401202 |
| | Buckles basic | | M | 401230 |
| | | | L | 401232 |
| | Magnet | | S | 401240 |
| | | | M | 401242 |
| | | | L | 401244 |
| | Ironclip | M | 401260 | |
| | | L | 401262 | |
| | Canvas | Magnet | M | 401400 |
| | | | L | 401402 |
| | | Ironclip | M | 401450 |
| L | | | 401452 | |

Colours



Black



Canvas

