

RESTRAINT BELTS THORAX BEDS 90



Description

Thorax belt for 90 cm bed.

Effects and benefits

The patient is immobilised in the abdominal, perineal and shoulder areas, preventing the patient from turning sideways and sitting up. It has been designed for optimal, comfortable and breathable distribution of the trunk area. Thanks to its design, maximum stability is guaranteed, allowing the user to move. Adaptable to any wheelchair.

Composition

This model is made of canvas material (75% cotton and 25% polyester), but can be made in polypropylene or quilted. The main difference in terms of the fabric is the washing temperature. The quilted fabric withstands up to 40°C and the polypropylene and canvas withstand up to 90°C.

As for the closure, the reference model has an ironclip closure.

Dimensions

- M: 85-109cm. This is the most popular size, but sizes are available:
- S: 60-90cm.
- L: 105-150cm.

Maintenance

- Maximum washing temperature 90°C.
- Do not dry clean.
- Do not tumble dry.
- Do not bleach.
- Do not iron.



Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Restraints

Restraint belts thorax beds 90

Models



Colours



Padded Navy
Blue



Canvas

| MODEL | CLOSURE | VERSION | FABRIC | REF. |
|--------------------------------------|----------|-----------|--------|--------|
| Restraint belts thorax beds 90 | Ironclip | Perineal | Padded | 402512 |
| | | | Canvas | 401414 |
| | | Abdominal | Padded | 402550 |

