

## RESTRAINT BELTS THORAX BEDS 90



### Description

Thorax belt for 90 cm bed.

### Effects and benefits

The patient is immobilised in the abdominal, perineal and shoulder areas, preventing the patient from turning sideways and sitting up. It has been designed for optimal, comfortable and breathable distribution of the trunk area. Thanks to its design, maximum stability is guaranteed, allowing the user to move. Adaptable to any wheelchair.

### Composition

This model is made of canvas material (75% cotton and 25% polyester), but can be made in polypropylene or quilted. The main difference in terms of the fabric is the washing temperature. The quilted fabric withstands up to 40°C and the polypropylene and canvas withstand up to 90°C.

As for the closure, the reference model has an ironclip closure.

### Dimensions

- M: 85-109cm. This is the most popular size, but sizes are available:
- S: 60-90cm.
- L: 105-150cm.

### Maintenance

Maximum washing temperature 90°C.

Do not dry clean.

Do not tumble dry.

Do not bleach.

Do not iron.



### Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

# Restraints

## Restraint belts thorax beds 90

### Models



### Colours



Padded Navy  
Blue



Canvas

MODEL	CLOSURE	VERSION	FABRIC	REF.
Restraint belts thorax beds 90	Ironclip	Perineal	Padded	402512
			Canvas	401414
		Abdominal	Padded	402550

