

POSITIONAL WEDGE LUMBAR



Description

Lumbar positional HR wedges for all types of wheelchairs.

Effects and advantages

Thanks to its shape, it prevents lateral trunk movements. Offers greater stability and support for the trunk. Dorsal-lumbar positioning orthosis that favors the correct positioning of the trunk, preventing the user from leaning sideways.

Composition

- Saniluxe outer cover: Fabric made up of two layers: 100% polyurethane lining and 100% polyester inner fabric. Fabric thickness 230g/m².
- Core: HR density 30 kg/m³.
- Closure system: polypropylene with Velcro.

Dimensiones

De la cuña: 33 cm de alto y 10 cm cada lado del triángulo

Maintenance

Simply wipe with a cloth or cloth. It can also be hand washed or machine washed at up to 90°C.

Tumble dry tolerant.


Do not bleach.



Warranty

The warranty period is two years from the date of purchase of the cushion. Products purchased on or after January 1, 2022 are covered by a three-year warranty. This warranty does not cover normal wear and tear of the product and does not replace. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in case of an incident in order to resolve it.

Lumbar positio Models



Colores



Graphite

MODEL	DIMENSIONS (CM)	REF.
Lumbar positional wedge	33 cm high and 10 cm each side of the triangle	705215

