TATAMI FLOOR AND SAFETY MATTRESS



Description

Tatami floor and safety mattress.

Effects and benefits

Cushions the impact in case of falling against the ground. Easy to store and takes up little space. Gentle unevenness on the sides of the mattress to facilitate the use of a wheelchair and facilitate the mobility of the user avoiding obstacles. It has rigid angles at the base to prolong the life of the mat and facilitate user access. The outer cover is very resistant, fire retardant, breathable, waterproof, stainresistant and antimicrobial.

Composition

Agglomerated tatami. Leatherette cover. PVC coating, applied over a 100% recycled polyester reinforcement, with a Nautic Plus finish.

Dimensions

Extended: 90x190x4cm Folded: 90x65x12cm

Maintenance

Hand or machine washable at 60°C with neutral detergent. It can be dry cleaned. Do not spin. Tumble dry at reduced temperature. Do not iron. Do not use bleach.






Warranty

The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the quality department of the company in case of incident to resolve it.

Geriatric

Tatami floor and safety mattress Models



Navy blue

MODEL	DIMENSIONS (CM)	COLOUR	REF.
Tatami floor and safety mattress	90x190x4	Navy blue	703025

