HAND REST





10 units of handrests for crutches.



Provides better grip and prevents the appearance of calluses.

Composition

Neoprene.

Dimensions

Ø6x14cm.

Maintenance

Wash periodically by hand or machine wash at a maximum temperature of 40°C.

Dry in a horizontal position without exposing to direct heat sources such as cookers, heaters, etc.

Warranty


The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.





Geriatric

Hand rest Models



Colours



MODEL	DIMENSIONS (CM)	REF.
Pack 10u hand rest	Ø6x14cm	801200

