

## SUAPEL SANITIZED SHEETS



### Description

Anti-decubitus sheets made of suapel sanitized. The polyester fibre exterior gives a feeling of comfort, being suitable for patients who need a feeling of warmth and extra comfort, thanks to its wool fleece fabric. It is suitable for use in wheelchairs, armchairs and beds.

### Consejos

The first prevention of bedsores in a sitting or lying position is to change the points of support of the body on the area where the body rests. On average, this should be done every three hours to reduce prolonged pressure on the bony prominences and to allow access of blood to the tissues.

The sheet should be checked to ensure that it fits the seat or bed where it is to be used.

### Instalación y precauciones

The shape of the suapel sheets is designed to adapt to the different types of seats and beds that patients may use and their needs.

Caustic products, paint strippers, solvents and sharp objects in direct contact with the product are prohibited. Avoid prolonged contact with liquids and moisture. Regularly inspect the sheet for normal wear and tear (changes in surface appearance, etc.).

### Composition

Exterior: polyester, terlenka.

### Dimensions

According to reference.

### Tratamiento del material y conservación

The sanitized cover treatment is applied to the surface of the fabric to create a highly effective anti-mite, anti-bacterial and anti-fungal barrier. A plus of protection for your bedding accessories.

Suapel sheets should preferably be stored flat and protected from natural light and humidity.

### Elimination

Do not burn or dispose of in nature or outside designated sites.

### Warranty

The warranty period is two years from the date of purchase of the cushion. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

Product lifetime: 1 year.

# Anti-decubitus

## Suapel sanitized sheets Models



MODEL	DIMENSIONS (CM)	REF.
Suapel sanitized sheets	40x60	105100
	50x70	105101
	60x80	105102
	70x80	105103
	60x90	105104
	70x90	105105
	80x90	105106
	90x150	105107
	110x170	105108
	70x190	105109
	80x190	105110
	90x190	105111
	105x190	105113
	120x190	105114
	135x190	105115
90x160	105117	
Suapel sanitized adjustable sheets	90x190x20	105125

### Colours



White



antiol or



anti bacterias



anti hongos



anti ácaros



anti CHINCHES



SUAVE



TRANSPIRABLE



LAVABLE