

## HEEL PAD BUCKLES (PAIR)



### Description

Leg restraints in different materials: padding, polypropylene and canvas with buckle, magnet or ironclip fastenings.

### Effects and benefits

Relaxed or firm fixation of one or both ankles.

### Composition

- Polypropylene: 100% polyester.
  - Padding: 20% suede (100% polyester), 20% velour (100% polyamide) and 60% foam.
  - Canvas: 75% cotton and 25% polyester.
- Buckle, magnet or ironclip fastening.  
100% polypropylene webbing.

### Dimensions

Fits any type of ankle.

### Maintenance

They do not require special care, just wash them periodically by hand or machine. Dry in a horizontal position without exposing to direct heat sources.

Maximum temperatures per material:

- Canvas and polypropylene: 90°C.
- Quilted: 40°C.

### Warranty


The warranty period is two years from the date of purchase. This warranty does not cover normal wear and tear of the product and does not replace legal warranties. Please contact the company's quality department in the event of an incident in order to resolve it.

# Restraints


## Heel pad buckles (pair)

### Models


01




02





03



04



05



MODEL	FABRIC	CLOSURE	REF.
Heel pad buckles	Polypropylene	Magnet	404300
		Buckles	404100
		Ironclip	404350
	Padded	Magnet	404500
		Buckles	404200
		Ironclip	404550
Canvas	Magnet	404400	
	Ironclip	404450	

### Colours



Black



Canvas

Polypropylene



Padded  
Navy Blue

